### **PRODUCT SPECIFICATION DC 16**

Edic.: 02



Register number: 01145-H **PRODUCT:** OPTIM-I.A.

**REVISION (yymmdd):** 

180530

#### **GENERAL CHARACTERISTICS:**

DESCRIPTION	Extender for boar semen in vitro at 15-17°C, for an approximate period of <b>4 – 5 days</b> depending on sperm quality that maintains its physic-chemical parameters (pH, pka, osmolality) as well as its spermatic functionality.  Product for Swine Artificial Insemination use.
PRESENTATION	
	Diluyente de media duración para semen de verraco Mid-term boar semen extender
PRODUCT PREPARATION FOR A CORRECT USE:	Dilute the content of the package at a proportion of <b>45,825 grams/litre</b> of bi-distilled or distilled water or, if failing, laboratory water of type II at a temperature of 37C°, stirring the solution to obtain a better result. Once the extender is reconstituted, mix with the ejaculate, trying that temperature differences between them be as low as possible.

#### **SPECIFICATIONS:**

PHARMACEUTICAL PRESENTATION FORM:	
	Physical form: odourless white powder.
	Characteristics of the packaging: heat-sealed aluminium complex in formats for 1 Litre, 5 Litres and 100 Litres of extender.
	Sealing system: heat-sealing.

**EP.SD.71** XX - Departamento vinculado (A: Almacén, SD: Diluyentes, TE: Taller, B: Bolsas)

YY - Número correlativo, unívoco e inequívoco



QUALITATIVE COMPOSITION:	
	Glucose.
	Sodium bicarbonate.
	Stabilizers of the spermatozoon plasmatic membrane.
	Excipients.
	It contains antibiotics complying with the European Directive 90/429/EEC.

MANUFACTURING METHOD:	
	MAGAPOR <b>OPTIM-I.A.</b> extender is manufactured with chemical products of pharmaceutical quality or, if failing, of food quality.
	There is no thermic process during the elaboration process, for a chemical synthesis of pure raw materials from the chemical or food industry is used, so there is no possibility of transmitting infectious agents. This product does not contain any protein of animal origin and thus is free of B.S.A. (bovine seroalbumin), preventing the transmission of potential pathogens.

OTHERS:	
METHOD OF ADMINISTRATION:	Not applicable (Do not inject. Do not ingest).
MAXIMUM PRESERVATION FOR USE AFTER RECONSTITUTION:	2-3 days - Store in a chamber between 2ºC y 8ºC, always checking that the product is hermetically sealed, protected from air and external agents.
POSSIBLE COLATERAL EFFECTS (LOCAL AND/OR GENERAL, INCOMPATIBILITIES AND PHARMACOLOGICAL ANTAGONISMS):	None.
INTOXICATION MAND OVERDOSE IN ANIMALS:	Not applicable.
INTOXICATION IN HUMAN BEINGS:	Not applicable (Do not ingest. Keep the product out of the reach of children).
UNDESIRED BIOLOGICAL EFFECTS:	Not applicable.
CONTROL OVER RESIDUES OF MEDICINAL PRODUCTS:	Not applicable.

# **PRODUCT SPECIFICATION**

## DC 16

Edic.: 02



CORRECT PRESERVATION OF THE PRODUCT:	Store in a cool and dry place, away from direct sunlight and moisture.
EXPIRY DATE (VALIDITY PERIOD):	
	24 MONTHS from manufacturing date (see package).

Page 3 of 3 **EP.SD.71**