


HYGITUBE

Disposable tube with filling needle for semen doses packaging devices.




Technical concept


Disposable tube with filling needle incorporated that avoids contamination problems in the packaging of semen doses.

» Depending on the needs of each user, **Hygitube can be discarded** when changing the boar, breed, working day, etc., with a maximum of 1,250 doses.

Advantages

- It avoids contamination.
- It guarantees that no remains of other ejaculates are packaged.
- It maximizes the dosing efficiency of the pump.
- It improves the insertion of the filling needle in Semenbag.
- It avoids sterilizing the silicone tubes.
- It ensures that all material in contact with the semen during the whole packaging process is disposable.

 20 cases
5 units/case

 15 Kg
40 x 60 x 54 cm

Characteristics

- Made of PVC and food grade silicone.
- Individually packaged.
- Easy installation.
- Compatible with Bagmatic, Simplymatic and Simplybag.
- It includes plastic needle.

